CLARIX EZ-4

PL

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Thank you for choosing the Tracer Clarix-EZ4. The device you have purchased is a breathalyser with a precise electrochemical sensor, used to measure the alcohol content in the exhaled air. Before using the device for the first time, please read the user manual below. This will enable you to take full advantage of the device's capabilities and use it correctly. We also recommend that you keep the manual in a place where it will be easily accessible if needed.

CAUTION! NEVER ALLOW YOURSELF OR ANOTHER PERSON TO DRIVE ANY VEHICLE UNDER THE INFLUENCE OF ALCOHOL. DOING SO ENDANGERS YOUR LIFE AND HEALTH AND MAY HAVE SERIOUS LEGAL CONSEQUENCES.

ADVANTAGES OF THE ELECTROCHEMICAL SENSOR

The electrochemical sensor detects alcohol molecules contained in exhaled air. The test devices in which it is used distinguish themselves from other types of breathalysers by their reliability and accuracy. Electrochemical sensors have a very fine-tuned sensitivity. The lifetime of such a sensor is also up to 5x longer than that of a solid-state sensor.


SAFETY INFORMATION FOR USE

Please read it in detail before you use it

1. Products can not be judged as the only device that gives guaranty whether it is allowed to drive after drinking alcohol.
2. Judgment is the general symptom and criteria. Even for the same alcohol concentration, intoxication will vary depending on physical condition and individual differences.
3. Regardless of whether the breathalyzer tester is used, the manufacturer and the sales dealer are not liable for drunk driving and legally prohibited behavior.
4. It takes about 20 minutes for alcohol to be absorbed from the digestive system into the blood, so it is recommended to measure alcohol for 20 minutes. If measured immediately after drinking alcohol, the measured results can reflect the alcohol concentration of the alcohol alcohol. Instead of the alcohol levels in the blood.
5. To obtain more accurate measurements, the interval between the secondary measurements must be greater than 3 minutes.
6. It must not be broken down or washed. Store it in places with less fragrance.
7. When measuring with low electric quality, there is a certain error between the measurement result and the actual value.
8. Do not store it in place on direct sunlight, do not make excessive collision or vibration. Stay away from the heating equipment.
9. Do not charge the breathalyzer near heat sources or on flammable surfaces.

10. Avoid excessive discharge – recharge the battery regularly to extend its lifespan.
11. Do not leave the device connected to the charger for a long time after it is fully charged.
12. Do not leave it in a heated car.

DEVICE DESCRIPTION

POWER ON/OFF:

1. Long press the middle key for 2 seconds, and there will be a „beep” for power-on
2. Default setting will shut down tester after 1 minute of no action. You can easily change it in the settings.


PREPARATION FOR TESTING:

1. Large number of experiments have shown that 15 to 20 minutes after drinking the alcohol, its content in the blood will reach the highest value, and then gradually decrease.
2. If tested immediately after drinking, the measurement result will increase significantly because there is still a small amount of high-concentration alcohol in the oral cavity.
3. Please take a deep breath before blowing to ensure that the continuous blowing time is greater than 6 seconds.
4. For the same amount of alcohol consumed, due to different personal constitutions, the measurement results may vary.

TESTING:

1. The screen shows „Start“, short press the middle key to start.
2. In the „Ready“ interface, short press the middle key to start.
3. After 10 seconds, it shows „Please Blow“, and a 30-second countdown will start.
 - When the countdown ends, it prompts „Timeout“ if the device will not notice proper blow.
4. The instrument monitors the blow air flow in real time, and the blow air time is about 6 seconds.
 - If the blow air time is less than 3 seconds, it prompts „Failure“.
 - If the blow air flow is too small, it prompts „Failure“.
 - If the blow air time is less than 6 seconds, it is marked as „Quick test“.
5. The long sound stops, stop blowing, and you may check results shown on the display.

- If the result is less than 20mg/100mL, it is on a green background.
- If the result is greater than 19mg/100mL, it is on a yellow background.
- If the result is greater than the alarm value, it is on a red background and continuous alarm with „beep, beep, beep”.

CHANGING UNITS:

After entering the settings mode, it is possible to set different units of measurement according to your requirements.

ALARM VALUE:

When the measurement result is greater than or equal to the set value, the buzzer emits a continuous „beep, beep, beep” sound and is displayed as a red background. This set value is defined as the alarm value.

ALARM VALUE SETTING:

When the screen shows „Alarm”, short press the middle key to enter the „Alarm value” mode. Short press the down key to decrease the alarm value, short press the middle key to increase the alarm value, and short press the back key to exit.

RECORD:

1. The pointer points to „Record”, and short-press the middle key to enter the „Record” interface.
2. Short-press the down key to scroll through the records.
3. Long-press the down key to scroll through the records quickly.
4. At the same time, long-press the return key to scroll through the records in the reverse direction.

BLUETOOTH PRINTING (OPTIONAL):

In the result interface, press the down key and then the middle key to print the current result.

In the record interface, press the middle key to print the current record.

MAINTENANCE CALIBRATION:

Do not place the instrument in extreme environments, such as long-term exposure to direct sunlight, high temperature, and high humidity.

Do not open the instrument housing. If it is opened, free maintenance services will not be provided.

According to the requirements of the international standard D/N31051, when the alcohol detector is used 500 times or for more than 6 months, it needs to be calibrated to ensure its accuracy.

Repair, maintenance, and calibration can only be provided by the company's service center!

IMPORTANT ADDITIONAL INFORMATION:

Some beverages, fruits, delicacies, and the surrounding air environment, etc. may affect the measurement results, so the measurement results are for reference only, and the seller and the manufacturer do not bear any legal responsibility to the user.

Our company strives to continuously improve the product performance, so our company reserves the right to modify and interpret this instruction manual.

The product is covered by the manufacturer's guarantee. The guarantee is only valid if the device is regularly serviced and calibrated at the manufacturer's service centre. At least once a year. The device is correctly calibrated and ready to use from new. Every breathalyser offered on the market requires periodic calibration, i.e. at least once every 12 months (we recommend every 6 months) the device must be configured to correctly indicate the amount of alcohol in the exhaled air. The need to calibrate the device is in no way a defect of the device and is not subject to warranty. If the device is not calibrated regularly it will most likely indicate distorted and incorrect measurement results.

Due to the technological sophistication of this type of sensor, a qualified service technician can readjust the sensor of the breathalyser so that it shows correct and reliable results.

Calibration is the process of readjusting the breathalyser so that it gives the most accurate results possible. The CLARIX EZ-4 breathalyser requires calibration every approx. 500 uses or every 12 months.

Calibration of the device within 24 months of purchase is free of charge. More information can be found at the address indicated: <https://www.tracer.pl/alkomat-tracer-ez-4>

Sensor	Elektrochemical
BAC measuring range	BAC: 0-500mg/100ml
Sensor warm-up time	20s
Operating conditions	-15°C / 40°C
Display	1.8" colour
Measurement memory	up to 1500
Battery	400mAh 3.7V built-in
Power supply	USB-C 5V 500mA
Dimensions	54*130*25mm
Net weight	100g
Measuring accuracy	BAC 0~50mg/100ml +-10%;
	BAC 50~100mg/100ml +-15%;
	BAC 100~200mg/100ml +-20%

Megabajt Sp. z o.o. hereby declares that the The Tracer Clarix EZ-4 radio breathalyser is, complies with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following website: www.tracer.pl/TRAALC47703

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